

## Recombinant bovine/porcine FGF2-G3 (145 aa) protein (Qk080-FG)



**Type:** Food grade proteins

**Available for purchase:** Unit Size (µg): 25, 50, 100, 500, 1000

Buy online with secure credit card or purchase order.

[View this product and buy online](#)

### Product Information

Recombinant bovine/porcine FGF2-G3 protein is a [thermostable engineered form of FGF-2](#) (bFGF) with an increased functional half-life from <10 h (wild-type) to >7 days (FGF2-G3). Bovine FGF2-G3 is used for the development of optimized serum-free culture media for species-specific bovine (cow) and porcine (pork) [cultivated meat](#) and veterinary research applications, including the expansion of bovine and porcine smooth muscle cells, bovine and porcine induced pluripotent, embryonic and mesenchymal stem cells.

Recombinant FGF2-G3 is used in B8 media ([Kuo et al. 2020](#)) and Beefy-9 media ([Stout et al. 2023](#)) for weekend free, high homogeneity induced pluripotent stem cell culture. FGF2-G3 also has applications in chemically defined stem cell and organoid culture media.

Qkine high quality food grade bovine/porcine FGF2-G3 has a molecular weight of 16.3 kDa. This protein is [animal origin-free](#), carrier-free and tag-free to ensure purity with exceptional lot-to-lot consistency.

#### Alternative protein names

Basic fibroblast growth factor, bFGF, FGF-β, FGF2, Fibroblast growth factor-basic, HBGF-2, FGF2-G3, FGF2-STAB, FGF2-STAB3

#### Molecular weight

16.3 kDa (monomer)

#### Protein Uniprot number

High purity thermostable bovine/porcine FGF2-G3

### **Species reactivity**

- bovine
- porcine

### **Product Information**

- High quality food grade recombinant protein
- >98%, by SDS-PAGE quantitative densitometry
- Expressed in *E. coli*
- Animal origin-free (AOF) and carrier protein-free
- Manufactured in the UK under a food manufacturing HACCP regime
- Lyophilized from Tris, NaCl, CyS, mannitol

### **Reconstitution instructions**

- Resuspend in sterile-filtered water at >50 µg/ml

### **Featured applications**

- Maintenance and expansion of bovine iPSC, ESC and primary cells
- Cellular agriculture and cultivated meat cell culture media optimization
- Cellular agriculture process development


### **Further quality assays**

- Mass spectrometry: single species with expected mass
- Recovery from stock vial: >95%
- Endotoxin: <0.05 EU/µg protein
- Full raw materials traceability, allergen analysis, CoO, CoA, beta-lactam-free and animal origin-free certification available

## Scientific Information


### Bioactivity

Recombinant bovine/porcine FGF2-G3 activity was determined using the Promega serum response element luciferase reporter assay in transfected HEK293T cells. EC<sub>50</sub> = 217 pg/ml (13 pM). Cells were treated in triplicate with a serial dilution of bovine/porcine FGF2-G3 for 3 hours. Firefly luciferase activity was measured and normalized to the control Renilla luciferase activity.



### Purity

Recombinant bovine/porcine FGF2-G3 145 aa migrates as a major band at approximately 16 kDa (monomer) in reduced (R) and non-reduced (NR) conditions. The dimeric form is also observed at approximately 32 kDa in the non-reduced condition. No contaminating protein bands are present. The purified recombinant protein (7 µg) was resolved using 15% w/v SDS-PAGE in reduced (+β-mercaptoethanol, R) and non-reduced (NR) conditions and stained with Coomassie Brilliant Blue R250.



**Original product page:** <https://ryan.calliope-alpha.ts.net/product/recombinant-bovine-porcine-fgf2-g3-145-aa-protein-qk080-fg/>

**PDF generated:** 12 May 2026

Copyright © 2026 by Qkine Ltd. All rights reserved including graphics and images.